

DATASHEET

TROLLEY TABLETOP

Product Number: 452219500831




Size:	2200*620*47 (L*W*T [mm])
X-Ray Attenuation:	< 1.1 Al. Eq. measured at 2.7 HVL
Material:	Sandwich of carbon fiber and Rohacell core, aluminum inserts

What makes our Trolley Tabletop unique?

- Lightweight (<23kg) tabletop for single side suspended patient tables
- Maximum patient load (line load) of 225kg
- Biocompatible surface
- Perfect X-ray homogeneity and low X-Ray attenuation
- 100% X-ray Product Test



Technical Drawing



NOTE: Standard color white, other colors available on request

Contact Information:

VDL Fibertech Industries
Hallenweg 15
5683 CT BEST
Netherlands
Phone: +31 (0)499 367676
Fax: +31 (0)499 367675